

TytoLogy2 functions


December 2015
Go Ashida (go.ashida@uni-oldenburg.de)

HPSearch2 - Main Window

The HPSearch2 Main Window is a complex software interface for neural data analysis. It features several panels and controls:

- Calibration:** Includes buttons for Load (L), Load (R), Plot CAL, and Delete. It shows file paths for L and R channels.
- TDT (Tactile Detection Threshold):** Contains a TDT Disable button, a dropdown for No_TDT, a Monitor OFF button, and a MonGain slider set to 1000.
- I/O Channels:** Controls for Input, Output (L), and Output (R) channels.
- TDT Settings:** Includes AcqDuration, SweepPeriod, TTLPulseDur, HP Freq, and LP Freq settings.
- Stimulus Settings:** Includes ISI, Duration, Delay, Ramp, RadVary, and FrozenStim options.
- Search:** A central panel with sliders for Latt, Ratt, ILD, ABI, BC%, Freq, BW, sAM%, and sAMHz. It includes a Search button and a checkbox for Left ON.
- Threshold Control:** Includes Threshold (Auto/Manual), Plot Y axis (Auto/Manual), Detection (Auto/Top/Bottom), and Sign (Plus/Minus) settings.
- Curves:** Includes a Stim type selector (BF, ITD, ILD, ABI, BC, FILD, Beat), a Side selector (Both, Left, Right), and a Run Curve button.
- Clicks:** Includes a Type selector (Cond(+), Rare(-), (+/-), (-/+)), a Side selector (Both, Left, Right), and a Run Clicks button.
- External Stimulus:** Includes a Side selector (Both, Left, Right) and a Run Stim button.
- Control:** Includes Abort, Pause, and Plot buttons.

HPSearch2 - Plot Window



Standard Flow of Experiments

* Before Starting Experiments

- + TDT hardware is already set-up and properly wired to the computer and other hardware.
- + Miniature microphones have been calibrated with **MicrophoneCal2**.
- + Microphone calibration data are saved as a Matlab file: **filename_fr2.mat**.

* Calibrating Earphones

- + Set the earphones (headphones) to the animal, open **HeadphoneCal2**, set the microphone calibration data file, set frequency ranges and other parameters, and calibrate the earphones.
- + Headphone calibration data are saved as a Matlab file: **filename_cal2.mat**.

* Starting HPSearch2c and TDT

- + Move to a directory where you want to save your data.
- + Open **HPSearch2c**.
- + Select your TDT hardware from the dropdown list on the **TDT** module of the GUI.
- + Check if your TDT hardware is turned on.
- + Start TDT by pressing the **[TDT Enable]** button.
- + Check if your TDT is running (see the processor usage on the your TDT hardware).

* Setting-up Calibration Files and Parameters

- + Load your earphone calibration file(s) using the **Calibration** module of the GUI.
- + Set I/O channels, TDT settings, Experiment Settings, Stimulus Settings in the GUI.
- + You may load your settings file to reuse your previous setting information.

* Running Search Stimuli

- + Turn on your earphones by clicking the **[Left ON]** and **[Right ON]** checkboxes in the **Search** module of the GUI.
- + Choose stimulus type and set frequency, intensity, ITD, etc.
- + Click the **[Search]** button to play sounds.
- + Sound stimuli are repeated until you stop them by clicking the **[Search]** button again.
- + Recorded neural data are shown in the plotting window, in which you can specify which panel to show.

* Setting Spike Discrimination Thresholds

- + Using the **Threshold Control** module, you can set the threshold of the built-in spike detector.
- + In the **Spike Analysis** module you can specify the analysis time range of your interest. Saved spike statistics is based on these parameters.

* Recording Neural Responses

- + The **Curve** module is for playing sounds (tone, noise, etc.) with varied parameters (frequency, ITD, intensity, etc.) and recording neural responses.
- + The **Click** module is for playing click stimuli and recording neural responses.
- + The **External Stimulus** module (beta version) is for playing external sound file and recording neural responses. You have to load your stimulus file before playing it.
- + Specify the desired parameters in these modules and **Stimulus Setting** module.
- + Click **[Run Curve]**, **[Run Clicks]**, or **[Run Stim]** button to start playing/recording.
- + Recorded data are saved in a Matlab file (only if the recording is not aborted) and in a binary file (**.dat2**).

* Closing HPSearch2c and TDT

- + By closing the main GUI window, the program automatically calls corresponding functions to stop TDT circuit and close the connection to the TDT hardware.

HPSearch2c Functions

* Some Important Notes

- + All parameters are stored in the structure named **handles**. To save parameters, you have to call **guidata(hObject, handles);** before returning from a sub-function.
- + Most TDT-related parameters are defined in **HPSearch2c_config** and stored under **handles.h2.config**.
- + Other parameters are defined in **HPSearch2c_init** and stored under corresponding entries of **handles.h2**.
- + In the following, RX6 is used as TDT hardware.
- + If other hardware is used, see corresponding functions defined in **HPSearch2c_config**.

* Initialization - Opening GUI & Setup

[Functions]	[Descriptions]
HPSearch2c.HPSearch2c_OpeningFcn	- This function is called when the GUI is opened.
HPSearch2c_Opening	- Setting up handles and parameters
HPSearch2c_config	- This function defines TDT-related parameters (e.g., circuit path and name, I/O functions)
HPSearch2c_init	- This function defines HPSearch2 parameters (e.g., search/curve parameters and limits, thresholds, etc.)

* Finalization - Closing TDT & Cleaning-up

[Functions]	[Descriptions]
HPSearch2c_CloseRequestFcn	- This function is called when the GUI is closed
HPSearch2c_Closing	- This function checks if TDT is running.
HPSearch2c_TDTclose	- If TDT is running, then the connection is closed.

* Abort Button

[Functions]	[Descriptions]
buttonAbort_Callback	- The on/off state of this button works as a flag for aborting.

- + The state of this button is monitored by Curve/Click functions.
- + If the button is detected to be ON, then the current recording is aborted.

* Pause Button

[Functions]	[Descriptions]
buttonPause_Callback	- The on/off state of this button works as a flag for stopping.

- + The state of this button is monitored by Curve/Click functions.
- + If the button is ON, then the current recording is stopped until the user presses the button again to make it OFF.

* Plot Button

[Functions]	[Descriptions]
buttonCallPlot_Callback	- This function calls an external plot function.
TytoView_simpleView	- Function for quickly seeing HPSearch data.

* Show Variable Button

[Functions]

buttonShowVal_Callback

[Descriptions]

- This function enables retrieving hidden variables.

* Search Button

[Functions]

buttonSearch_Callback

HPSearch2c_enableUIs

HPSearch2c_Search

***** Setting TDT parameters (duration, delay, etc.) *****

HPSearch2c_RX6settings

HPSearch2c_RX8settings

RPsamplefreq

(Location = C:\TytoLogy\toolbox\TDT\Functions\RP\RPsamplefreq.m)

RPsettag

(Location = C:\TytoLogy\toolbox\TDT\Functions\RP\RPsettag.m)

***** Making stimulus waveforms and setting the attenuators *****

HPSearch2c_searchParamFromUI

TytoLogy2_mergecal

TytoLogy2_interpcal

syn_headphone_tone

syn_headphone_noise

syn_headphone_amtone

syn_headphone_amnoise

computeLRspl

TytoLogy2_figureAtten

PA5setatten

***** Playing sounds and recording data *****

HPSearch2c_spikeio

RPwriteV

RPtrig

RPfastgettag

RPgettag

RPreadV

*** Plotting data ***

HPSearch2c_plotParamFromUI

HPSearch2c_spikedetect

HPSearch2c_plotResponse

[Descriptions]

- This function is called when the [\[Search\]](#) button is pressed.

- Enabling/disabling GUI

- Main function for playing search stimuli

- Defined as **handles.h2.config.TDTsetFunc.**

- Internally called by **HPSearch2c_RX6settings.**

- This function gets the sampling frequency of the TDT hardware.

- This function sets various parameters to TDT hardware.

- Getting search parameter values from the GUI.

- Merging L and R calibration data.

- Interpolating calibration data to fit the stimulus frequency profile

- Making tonal stimuli

- Making noise stimuli

- Making AM tone stimuli

- Making AM noise stimuli

- Calculating sound intensities (SPL) from ABI and ILD

- Calculating attenuation levels from SPL and calibration data

- Defined as **handles.h2.config.setattenFunc.**

Setting attenuation to TDT PA5 hardware

- Defined as **handles.h2.config.ioFunc.**

Playing sounds and recording neural responses.

- Setting output stimulus to TDT

- Sending trigger command to TDT

- Checking if stimulus sweep ends

- Getting index of recorded data from TDT

- Getting recorded data from TDT

- Getting plot/threshold settings from the GUI

- Spike detection using the threshold settings

- Plotting traces, rasters, PSTH, ISIH, etc.

* Curves Button

[Functions]	[Descriptions]
buttonCurve_Callback	- This function is called when [Run Curve] button is pressed.
HPSearch2c_enableUIs	- Enabling/disabling GUI
HPSearch2c_Search	- Main function for recording
***** Setting TDT parameters (duration, delay, etc.) *****	
HPSearch2c_RX6settings	- Defined as handles.h2.config.TDTsetFunc .
HPSearch2c_RX8settings	- Internally called by HPSearch2c_RX6settings .
RPsamplefreq	- This function gets the sampling frequency of the TDT hardware.
RPsettag	- This function sets various parameters to TDT hardware.
***** Setting calibration data and checking frequency and intensity parameters *****	
TytoLogy2_mergecal	- Merging L and R calibration data.
TytoLogy2_interpcal	- Interpolating calibration data to fit the stimulus frequency profile
***** Setting Loop variables (freq, ITD, etc.) *****	
-- Data are stored in loopvars .	
***** Setting output file names *****	
TytoLogy2_buildFileName	
***** Storing stimulus parameters and randomizing order *****	
-- Data are stored in stimcache .	
***** Making stimulus waveforms *****	
syn_headphone_tone	- Making tonal stimuli
syn_headphone_noise	- Making noise stimuli
syn_headphone_amtone	- Making AM tone stimuli
syn_headphone_amnoise	- Making AM noise stimuli
***** Calculating attenuation levels *****	
computeLRspl	- Calculating sound intensities (SPL) from ABI and ILD
TytoLogy2_figureAtten	- Calculating attenuation levels from SPL and calibration data
***** Setting up binary data file *****	
TytoLogy2_writebinary	- Writing headers of binary data file
***** Getting spontaneous response data *****	
PA5setatten	- Defined as handles.h2.config.setattenFunc .
HPSearch2c_spikeio	- Defined as handles.h2.config.ioFunc .
***** Main loop: playing sounds and recording data *****	
PA5setatten	- Defined as handles.h2.config.setattenFunc .
HPSearch2c_spikeio	- Defined as handles.h2.config.ioFunc . Playing sounds and recording neural responses.
RPwriteV	- Setting output stimulus to TDT
RPtrig	- Sending trigger command to TDT
RPfastgettag	- Checking if stimulus sweep ends
RPgettag	- Getting index of recorded data from TDT
RPreadv	- Getting recorded data from TDT
TytoLogy2_writebinary	- Writing binary output data
*** Plotting data ***	
HPSearch2c_plotParamFromUI	- Getting plot/threshold settings from the GUI
HPSearch2c_spikedetect	- Spike detection using the threshold settings
HPSearch2c_plotResponse	- Plotting traces, rasters, PSTH, ISIH, etc.
HPSearch2c_plotCurve	- Plotting tuning curve
*** Check pause and abort flags, stop to fit ISI, and repeat ***	
-- If the abort flag is set, then recording is aborted.	
*** Saving and plotting ***	
TytoView_simpleplot	- Plotting recorded data (for quick checking)

* Clicks Button

[Functions]	[Descriptions]
buttonClick_Callback	- This function is called when [Run Clicks] button is pressed.
HPSearch2c_enableUIs	- Enabling/disabling GUI
HPSearch2c_Click	- Main function for playing clicks and recording
***** Setting TDT parameters (duration, delay, etc.) *****	
HPSearch2c_RX6settings	- Defined as handles.h2.config.TDTsetFunc .
HPSearch2c_RX8settings	- Internally called by HPSearch2c_RX6settings .
RPsamplefreq	- This function gets the sampling frequency of the TDT hardware.
RPsettag	- This function sets various parameters to TDT hardware.
***** Setting output file names *****	
TytoLog2_buildFileName	
***** Setting loop variables (reps, ITD, etc.) *****	
-- Data are stored in loopvars .	
***** Setting up binary data file *****	
TytoLog2_writebinary	- Writing headers of binary data file
***** Getting spontaneous response data *****	
PA5setatten	- Defined as handles.h2.config.setattenFunc .
HPSearch2c_spikeio	- Defined as handles.h2.config.ioFunc .
***** Main loop: playing sounds and recording data *****	
syn_headphone_click	- Making tonal stimuli
PA5setatten	- Defined as handles.h2.config.setattenFunc .
HPSearch2c_spikeio	- Defined as handles.h2.config.ioFunc .
	Playing sounds and recording neural responses.
RPwriteV	- Setting output stimulus to TDT
RPtrig	- Sending trigger command to TDT
RPfastgettag	- Checking if stimulus sweep ends
RPgettag	- Getting index of recorded data from TDT
RPreadv	- Getting recorded data from TDT
TytoLog2_writebinary	- Writing binary output data
*** Plotting data ***	
HPSearch2c_plotParamFromUI	- Getting plot/threshold settings from the GUI
HPSearch2c_spikedetect	- Spike detection using the threshold settings
HPSearch2c_plotResponse	- Plotting traces, rasters, PSTH, ISIH, etc.
HPSearch2c_plotCurve	- Plotting tuning curve
*** Check pause and abort flags, stop to fit ISI, and repeat ***	
-- If the abort flag is set, then recording is aborted.	
*** Saving and plotting ***	
TytoView_clickplot	- Plotting recorded click response data (for quick checking)

* Run Stim Button (External Stimulus module)

[Functions]	[Descriptions]
buttonExtStimRun_Callback	- This function is called when [Run Stim] button is pressed.
HPSearch2c_enableUIs	- Enabling/disabling GUI
HPSearch2c_RunStim	- Main function for playing external stimulus and recording
***** Setting output file names *****	
TytoLogy2_buildFileName	
***** Setting TDT parameters (duration, delay, etc.) *****	
HPSearch2c_RX6settings	- Defined as handles.h2.config.TDTsetFunc .
HPSearch2c_RX8settings	- Internally called by HPSearch2c_RX6settings .
RPsamplefreq	- This function gets the sampling frequency of the TDT hardware.
RPsettag	- This function sets various parameters to TDT hardware.
***** Setting calibration data and checking frequency and intensity parameters *****	
TytoLogy2_mergecal	- Merging L and R calibration data.
***** Setting up binary data file *****	
TytoLogy2_writebinary	- Writing headers of binary data file
***** Setting attenuator *****	
PA5setatten	- Defined as handles.h2.config.setattenFunc .
***** Getting spontaneous response data *****	
HPSearch2c_spikeio	- Defined as handles.h2.config.ioFunc .
***** Main loop: playing sounds and recording data *****	
PA5setatten	- Defined as handles.h2.config.setattenFunc .
HPSearch2c_spikeio	- Defined as handles.h2.config.ioFunc . Playing sounds and recording neural responses.
RPwriteV	- Setting output stimulus to TDT
RPtrig	- Sending trigger command to TDT
RPfastgettag	- Checking if stimulus sweep ends
RPgettag	- Getting index of recorded data from TDT
RPreadv	- Getting recorded data from TDT
TytoLogy2_writebinary	- Writing binary output data
*** Plotting data ***	
HPSearch2c_plotParamFromUI	- Getting plot/threshold settings from the GUI
HPSearch2c_spikedetect	- Spike detection using the threshold settings
HPSearch2c_plotResponse	- Plotting traces, rasters, PSTH, ISIH, etc.