

How to Install TytoLogic and TytoLogic2

Go Ashida (ashida@umd.edu)
June 2012

How to Install

- (1) Copy the **TytoLogic** folder to your hard drive: **C:\TytoLogic**.
If you already have TytoLogic, you can skip this step.
- (2) Copy the **TytoLogic2** folder to your hard drive: **C:\TytoLogic2**.
- (3) Start your Matlab and call **setrigpath2** to set your PATH to TytoLogic and TytoLogic2.
Now you are ready to use TytoLogic2.

Notes

- * Since TytoLogic2 internally uses several functions of TytoLogic, you should install both TytoLogic and TytoLogic2 even if you use only TytoLogic2.
- * You may install TytoLogic/TytoLogic2 to any folder other than **C:\TytoLogic** or **C:\TytoLogic2**, but you will have to edit the following Matlab script files so that the program can find corresponding scripts and circuits properly.
 - setrigpath2.m
 - MicrophoneCal2_init.m
 - HeadphoneCal2_init.m
 - HPSearch2_config.m
- * For technical reasons, TytoLogic and TytoLogic2 do not use the same data structure. You may need to re-calibrate your microphones and earphones/headphones using MicrophoneCal2 and HeadphoneCal2, before using HPSearch2.